

# CHLAMYDIA CEL Pn

#### INTENDED USE AND PRINCIPLE OF THE TEST

The Chlamydia Cel Pn If Test is an *in vitro* indirect immunofluorescence test for the direct detection of *Chlamydia pneumoniae* organisms in clinical specimens. A monoclonal antibody to *C. pneumoniae* binds specifically to the organisms present in fixed specimens and a second antibody, FITC-conjugated goat antimouse lg, stains the *C. pneumoniae* organisms.

#### CONTENTS OF THE KIT

RC RN

R٨

The Chlamydia Cel Pn range is available in a number of formats, including: 10 Test Demo (KC3 Demo), 50 Test (KC3), 200 Test Bulk, Positive Control Slide only, and Mounting Fluid only. Note that the Chlamydia Cel Pn MAb Reagent and Anti-mouse Ig-FITC Reagent has been optimised for use with Cellabs *Chlamydia pneumoniae* Positive Control Slide and Mounting Fluid.

		KC3 Demo	KC3 Standard	Bulk
C3	Chlamydia Cel Pn MAb Reagent	0.25mL	1.25mL	5mL
Μ	Anti-mouse Ig-FITC Reagent	0.25mL	1.25mL	5mL
	Positive Control Slide (Single Use only)	1	1	-
MF	Mounting Fluid	2.5mL	2.5mL	-
	Tests	10	50	200


All components should be stored at 2-8<sup>o</sup>C, and are supplied ready for use. Expiry dates are clearly marked on each kit component and on the box and do not change once opened.

### MATERIALS REQUIRED BUT NOT PROVIDED

- Microscope slides with 6-8mm diameter wells
- Acetone or methanol or acetone for specimen fixation
- Precision pipette for delivering 25µL
- Humid chamber
- Wash bath
- Phosphate buffered saline (PBS) for washing step
- Cover slips
- Non-fluorescing immersion oil
- Fluorescence microscope with filter system for FITC (maximum excitation wavelength 490nm, mean emission wavelength 530nm) and x600-x1000 magnification

#### **DIAGRAM FOR USE**

Use Cellabs Instructions for Use Insert contained in kit when performing test, and refer to Material Safety Data Sheet (MSDS) for further information.





# **Cellabs Product Profile**

#### READING AND INTERPRETATION OF RESULTS AND DIAGNOSIS

*C. pneumoniae* elementary bodies and reticulate bodies fluoresce bright apple green (with typical morphology) against a background of reddish/brown counterstained material. Results should be compared with the Positive Control Slide.

#### PERFORMANCE DATA FOR CHLAMYDIA CEL Pn

#### Sensitivity/Specificity

А	n = 53. Chlamydia Cel Pn versus culture.	Sensitivity: 100%		
		Specificity: 100%		

#### **Repeatability & Reproducibility**

5 identical antigen samples were tested by 3 different operators and there was 100% correlation in the results.

**Cross reactivity** The Chlamydia Cel Pn does not cross react with:

Staphylococcus epidermidis Staphylococcus aureus Beta haemolytic Streptococcus group A Streptococcus pneumoniae Haemophilus influenzae Haemophilus parainfluenzae Escherichia coli Klebsiella pneumoniae Proteus mirabilis Pseudomonas aeruginosa Bacteroides fragilis Viridans Streptococcus Candida albicans Mycoplasma pneumoniae Chlamydia trachomatis Chlamydia psittaci Respiratory Syncytial Virus Adenovirus Parainfluenzae type 3 Influenzae type B Cytomegalovirus

## For Ordering Assistance:

See Your Local Distributor:

OR

Cellabs Pty Ltd Unit 7, 27 Dale Street (PO Box 421) Brookvale, NSW 2100 Australia Tel: +61 2 9905 0133 Fax: +61 2 9905 6426 Web: http://www.cellabs.com.au Email: sales@cellabs.com.au



Authorised Representative in the European Community: MDCI Ltd Arundel House 1 Liverpool Gardens Worthing, West Sussex BN11 1SL United Kingdom