SELECTROL® – ENGLISH

Selectrol® discs are intended for use in microbiological laboratories for the control of test methods. Being **1st generation derivatives** traceable to recognised culture collection strains, Selectrol® discs are acceptable in accredited laboratories for the production of **working stock** cultures (see Limitations).

Selectrol® discs are freeze-dried microbial preparations manufactured exclusively from NCTC® (National Collection of Type Cultures) and NCPF® (National Collection of Pathogenic Fungi) cultures. The viability of the organisms is stabilised by the unique formulation of the disc.

Precautions & storage

- Selectrol® discs contain viable pathogenic organisms and must only be used in appropriately equipped laboratories by microbiologists or persons under the supervision of microbiologists qualified by training or experience to work with micro-organisms.
- Store at the temperature indicated on the label. Loss of viability may occur if the discs are left above this temperature range for longer than necessary.
- Allow the vial to reach room temperature before opening and replace the cap on each vial immediately after a disc has been removed. Moisture in air entering the vial will cause a reduction in the number of viable organisms, eventually leading to complete loss of viability.
- Do not use discs that are past their expiry date. These may show a loss of viability.

Removing the disc from the vial: A disc is easily removed from the vial by using sterile forceps or a sterile 10µl inoculating loop.

For use on solid media: Place a disc on appropriate solid medium. Allow the disc to soften for 10–15 minutes. The plate may be placed in an incubator to accelerate the process. Spread the softened disc around the plate and incubate under optimum conditions for the strain.

For use in liquid media: Place the disc in 1-10mls of the appropriate broth. Avoiding the production of aerosols, shake the broth gently to dissolve the disc and incubate under optimum conditions. Alternatively, for rapid use, dissolve the disc, allow the culture to recover by incubating at 35–37°C for 1 hour and then use immediately. By experimentation you will be able to find the best dilution for your own application.

Limitations: Repeated subculture can cause the characteristics of a strain to change. A Selectrol® disc is a first generation subculture from a **master culture** sourced from Public Health England Culture Collections, and is designed to be used to obtain **working stock** cultures for use in testing. It is generally accepted that no more than five passages (successive subcultures) should be made from the **master culture**, in order to avoid genetic drift and mutant selection. Therefore, no more than four passages from the Selectrol® **working stock** culture should be made.

Breakages, spillage & disposal

- All exposed discs, contaminated packaging and broken glassware should be placed in a suitable container and either incinerated or autoclaved at 121°C for 30 minutes.
- Unbroken vials may be removed with forceps, washed in a suitable bactericidal solution, rinsed, dried and retained for further use.
- All contaminated surfaces must be disinfected with a suitable bactericidal solution.

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Using the On-line Selectrol® Certificate System

Access the TCS Biosciences Ltd website by entering: www.tcsbiosciences.co.uk into the address box of your browser. Why not add us to your favourites?

The Selectrol® certificates of analysis button, located in our navigation bar, will take you directly to the On-line Selectrol® Certificates System.

Once the main page has loaded you may search for certificates using the following criteria:

- Code
- Lot Number
- Strain
- Strain Designation

SEARCH TIPS

Searches can use % as a 'wildcard' character, e.g. 12%0.

If you would like to retrieve more than one certificate, multiple search criteria and combinations can be used, provided that entries are separated by a comma. For example: 09800101, 09800102, or Rhodococcus, Streptococcus, Aspergillus

Database listings can be sorted by selecting the desired sort option from the drop down menu located to the right of the search box.

The system automatically displays 250 results per page, this number can be changed by selecting an option from the drop down menu at the bottom of the page and clicking 'Items per Page' button.

Once your search has been successful, pages can be advanced or reversed by clicking one of the links at the bottom of the list, e.g.:

- | | |
|---------------|--|
| [1] | Takes to you to the first page of the list |
| [2], [3], ... | Takes you to subsequent pages |
| [Next] | Advances to the next page of the list |
| [Last] | Takes you to the last page of the list |
| [Show All] | Displays all results. |

VIEWING CERTIFICATES

Clicking 'Display' opens the certificate of choice.

Clicking 'Print' will open the print dialogue box on your computer, allowing you to print the certificate to your chosen printer.

Clicking 'Back to List' will return you to the search results screen.

CLEARING THE STORED SEARCH TERMS


Click 'Reset' to reset the search terms so that all the records are available for searching or viewing (i.e. starting all over again).

t: +44 (0)1296 714222 f: +44 (0)1296 714806 e: sales@tcsgroup.co.uk i: www.tcsbiosciences.co.uk

Selectrol® is endorsed by the Public Health England (PHE) Culture Collections

PHE Culture Collections and TCS Biosciences Ltd have agreed a manufacturing licence with the result that Selectrol® will continue to be manufactured exclusively from recognised source strains from the PHE Culture Collections.

This collaboration reinforces the quality of the Selectrol® range of micro-organisms while ensuring all the benefits of first generation cultures.

Selectrol® discs are manufactured exclusively from NCTC® (National Collection of Type Cultures) and NCPF® (National Collection of Pathogenic Fungi) cultures and are guaranteed first generation derivatives from the source strain.

ATCC® strains are listed as a reference only. ATCC® is a registered trademark of the American Type Culture Collection.

Selectrol® Quality Control Testing

Selectrol® batches are tested for a range of identification and characterisation attributes in our UKAS accredited testing laboratory number 2496. A test report for each batch of Selectrol® can be accessed via our website, www.tcsbiosciences.co.uk. The reporting of the test results via the website forms part of our UKAS accreditation.

Please contact our Customer Service team for a copy of our UKAS certificate and schedule of accreditation. Our Schedule can also be downloaded from the UKAS website – www.ukas.com

Selectrol®-skiver er tiltenkt bruk i mikrobiologiske laboratorier for kontroll av testmetoder. Da Selectrol®-skivene er **1. generasjons derivativer**, sporbare til hetteglass med kontrollstammer, kan de aksepteres i akkrediterte laboratorier for fremstilling av kontrollstammekulturer (se Begrensninger).

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