

Certificate of Analysis

Product Description: Donor Horse Serum

Country of Origin: United Kingdom

Batch No: 22981

Expiry: 2015.08

Sterility Testing

Bacteria	Not Detected
Yeasts and Fungi	Not Detected
Mycoplasma	Not Detected
Viruses:	
Equine Infectious Anaemia	Not Detected
Equine Viral Arteritis	Not Detected

Physical and Biochemical Analysis

pH	7.59
Osmolality	288 mosmol/kg
Haemoglobin	0.169 mg/ml
Total Protein	67.9 mg/ml
Endotoxin	0.272 EU/ml
Visual check	Satisfactory

Functional Testing

Cell Growth (L929)	103% of control
Cell Growth (SP2)	99% of control
Plating Efficiency (L929)	95% of control

No evidence of cytotoxicity noted

Documentation approval

Quality Assurance – Date

Product Release

Quality Control - Date

TCS Biosciences Ltd

Botolph Claydon, Buckingham MK18 2LR

t +44 (0)1296 714222

f +44 (0)1296 714806

e sales@tcsgroup.co.uk

i www.tcsbiosciences.co.uk