

Product Note

Cryptosporidium & Giardia

- EasySeed™ is a quality control reference sample containing 100 *Cryptosporidium* and 100 *Giardia* (C&G) per tube
- Each kit contains 10 EasySeed™ tubes
- *Cryptosporidium* and *Giardia* in EasySeed™ are **inactivated** – safe to handle
- 4 month shelf life

Why use EasySeed™?

- Precise number of C&G allows accurate measurement of the test recovery
- Easiest way to meet **regulatory** Quality Control procedures
- **Saves labor and time** in preparing reference samples
- Inactivated organisms make it **safe** to handle in the laboratory
- **Longer shelf life** than live organisms allows flexibility in work schedules

Approvals

- United States Environmental Protection Agency
- United Kingdom Drinking Water Inspectorate
- Japan Water Works Association