

EasySeed™

- EasySeed™ is a quality control reference sample containing 100 *Cryptosporidium* and 100 *Giardia* (C&G) per tube
- Each kit contains 10 EasySeed™ tubes
- *Cryptosporidium* and *Giardia* in EasySeed™ are **inactivated** – safe to handle
- 4 month shelf life

Why use EasySeed™

- **Precise** number of C&G allows accurate measurement of the test recovery
- An easy way to meet **regulatory** Quality Control procedures
- **Saves labor and time** in preparing reference samples
- Inactivated organisms make it **safe** to handle in the laboratory
- **Longer shelf life** than live organisms allows flexibility in work schedules

Approvals

- United States Environmental Protection Agency
- United Kingdom Drinking Water Inspectorate
- Japan Water Works Association



EasySeed™ packaging



Distributed by:
TCS Biosciences Ltd
Botolph Claydon, Buckingham, MK18 2LR, United Kingdom
t: +44 (0) 1296 714222, f: +44 (0) 1296 714806, e: sales@tcsgroup.co.uk, i: www.tcsbiosciences.co.uk

